












Sr No	Area	Model	Product Details	Picture
1	CAUTERY	SMART 4 Model 306	Smart 4 is a micro processor controlled based model for surgeries such as Electro vaporization Gastroenterology, cardiology, Proctology, Gynecology, Urology, Neurology, Orthopaedics & many other type of surgeries.	
2	CAUTERY	400 W DIGITAL CAUTERY MPM 333	Tanvi MPM 333 is micro controller based model for surgeries such as Electro vapourization Gastroenterology, Cardiology, Proctology, Gynaecology, Urology, Neurology, Orthopaedics & many other types of surgeries. Digital Power adjustment for separate modes Feather touch keys are used for power setting levels in all operating mode Audio visual Alarm for disconnecting of patient plate Tissue response circuit is provided Programme selector switch is provided	
3	CAUTERY	300 W DIGITAL CAUTERY MPM 3033	We offer to our clients "Tanvi 330" units which are designed to perform both simple procedures and highly complex surgeries. Further they also provide a smooth cut, minimum tissue damage and effective coagulation.	
4	CAUTERY	400 W ANALOG CAUTERY MPM 303	Solid State surgical Diathermy 400 w with Handswitch Pencil Analogue type	
5	CAUTERY	300 W ANALOG CAUTERY 300 B	Solid State surgical Diathermy 300 w with Handswitch Pencil Analogue type	
6	CAUTERY	Radio Frequency Cautery (3MHz)	We offer to our clients superior range of Radio Frequency Cautery- 3 , which delivers radio frequency (3MHz) wave for cutting & coagulation of soft tissues. These radio waves are used for vaporizing the cells without applying any pressure on the tissue.	
7	OT Light	Alma Series - Tanvi 304	4 reflector (heart shape)	
8	OT Light	Alma Series - Tanvi 3044	4 + 4 reflector (heart shape)	
9	OT Light	Alma Series - Tanvi 305	4 reflector (Pentagonal shape)	
10	OT Light	Alma Series - Tanvi 3055	4 + 4 reflector (Pentagonal shape)	
11	OT TABLE	C-ARM COMPATIBLE Hydraulic Table (3015)	Tanvi c-arm semi motorized table is designed in accordance with orthopedic surgeons of radiographic capability & general surgery. Up and Down Movement is controlled by feather touch buttons on remote control without disturbing operating team. Radiolucent top allows clear image on c-arm machine	
12		C-ARM COMPATIBLE Semi Motorized OT table (3015M)	Tanvi c-arm semi motorized table is designed in accordance with orthopedic surgeons of radiographic capability & general surgery. Up and Down Movement is controlled by feather touch buttons on remote control without disturbing operating team. Radiolucent top allows clear image on c-arm machine monitor Eccentric table top allows max coverage for c-arm movement head & leg portion is reversible Floor locking is achieved by single knob or by a single pedal system Except up and down movement all other movements are achieved by rotating individual gear	