

Sr No	Area	Model	Product Details	Picture
1	Slides	Obio	A standard microscope slide measures about 75 mm by 25 mm (3" by 1") and is about 1 mm thick. A range of other sizes are available for various special purposes, such as 75 x 50 mm for geological use, 46 x 27 mm for petrographic studies, and 48 x 28 mm for thin sections.	
2	Micro Tips	Tarsons	Medical grade polypropylene autoclavable, free of lubricants, dyes and heavy metals. Certified DNase, RNase, Human DNA and PCR inhibitor free. Sterile Tips are Pyrogen Free.	
3	Vacutainer	BD	Tubes contain spray-coated silica and a polymer gel for serum separation. They are used for serum determinations in chemistry. BD Vacutainer® SST™ Tubes provide an efficient means of serum sample preparation and help to improve laboratory workflow.	
4	Urine Container	Obio	It is typically used. A special pan that fits in the toilet or a urinal may be used to collect the urine. You will need to transfer the urine from the collecting container to the storage container. You will need to keep it cold.	
5	Double Cap Vials	Obio	It is basically used for draw the blood from the patients. It is basically 3 types EDTA, Clot and Sugar.	
6	Test Tubes	Obio	A test tube, also known as a culture tube or sample tube, is a common piece of laboratory glassware consisting of a finger-like length of glass or clear plastic tubing, open at the top and closed at the bottom.	
7	Pipette Stand	Obio	Pipettes stand is basically used for hold the multiple pipettes in. the single stand	



OBio...

Supporting Care to Cure

A Division of Orchard Medical